

AK-C

Liner de Cerámica



AK
MEDICAL

CE 2292



CARACTERÍSTICAS

Material Cerámica Delta.

Compatible con copas AK-AC-II-TP, AK-AC-II-TTM-I, AK-AC-II-TTM-V.

AK-C Liner de Cerámica

DESCRIPCIÓN	REF	CÓDIGO	DIÁMETRO (mm)	COMPATIBILIDAD CON CABEZA FEMORAL
AK-C Liner de Ceramica	28/36	A2400-2836	36	28
	28/38	A2400-2838	38	
AK-C-Liner Ceramic	32/40	A2400-3240	40	32
	32/42	A2400-3242	42	
Material: Ceramic Delta	36/44	A2400-3644	44	36
	36/46	A2400-3646	46	
	36/48	A2400-3648	48	
	36/50	A2400-3650	50	
	36/52	A2400-3652	52	
	36/54	A2400-3654	54	